

Quality of life matters

CACVUS MEDICAL
Stores Limited

TruHemo™ POCT Coagulation Analyzer



 **One Drop of
Whole Blood**



**PT/INR APTT
FIB TT ACT**

Portable

Fast

Accurate



CE FSC



MHRA

IVD



Advantage

- Plug and play, easy to use
- Whole blood test
- Accuracy and Stable
- Fast results within 3 mins
- Room temperature for 18 months



How to use



1 Fingerstick microsample collection.



2 Add the sample to the reagent card.

Specification

Storage Conditions	Ambient temperature: -10°C to 50°C Relative humidity : ≤80%	Service Life	5 Years
Power Supply	AC220V	Test Conditions	Detection Temperature: 15°C to 35°C
Weight	Approx. 1.0 kg		Humidity Range: 10%~80% RH
Dimensions	233mm(L) × 197mm(W) × 94 mm(H)		Atmospheric Pressure Range: 70~106 kPa

Performance parameters

Analyte	PT	APTT	FIB	TT	ACT
Measurement Range	7~90S	18~120S	1.5~6g/L	10~90S	50~500S
Pre heating time	<40s				
Sample Volume	20 uL				
Sample Type	Whole Blood				



P.O BOX 73360, KAMPALA
Liberty Towers, Plot No. 59A-59B,
level 4, Room No.L4A-01. K'Ila Road.
Tel: 0393194409 / 0709760701
Email: cactusmedicalstoresltd@gmail.com

