

MK-A06

Hospital Accompanying Sleeping Chair



*Better care
Better comfort!*

MK-A06

HOSPITAL ACCOMPANYING SLEEPING CHAIR

DESCRIPTION

This accompanying chair has a simple design and can be converted into a simple bed. It has wheels for easy transportation and the overall material is easy to clean, making it an indispensable configuration in the ward.

- FEATURES
- Ergonomic design,with armrest and headrest.
 - Can be unfolded into a bed and retracted into a chair.
 - One-handed operation, convenient and fast.
 - Designed with wheels for easy transportation.

STANDARD SPECIFICATION

Model	MK-A06
Folded size	D670*W590*H910mm
Unfolded size	L1830*W590*H600mm
Castor diameter	Dia 40mm



Material and structure

Steel with epoxy painted frame, foam mattress covered with leather, easy to clean.



Design

It can be pulled open to become a companion bed, making it convenient to accompany the patient.



Castor

There are 6 castors with rubber feet, which are convenient for deployment, retrieval and transportation, while ensuring stability.



Pull cord design

The pull cord design allows you to retract the accompanying chair with one hand.



Humanized design

With armrests and headrest, convenient for caregivers and improved safety.