

YA-D5-10

Electric Hospital Bed



*Better care
Better comfort!*

YA-D5-10

Electric Hospital Bed

DESCRIPTION

The five-function electric hospital bed with soft connection, and the adjustment part has a soft connection to improve the comfort of the patient. Central lock, the detachable headboard and footboard of the bed greatly provide convenience for medical personnel.

- FEATURES
- Four split side rail, drop down under the bed platform
 - Hand remote control
 - Detachable head/footboard with easy lock
 - With articulated ABS soft joint between different platform sections
 - Central brake system with a wheel diameter of 125mm
 - Auto-counter to provide the patient from sliding onto foot end
 - Manual CPR on both sides(optional)

STANDARD SPECIFICATION

Overall size	2160*1050 *450-700mm
Lying size	1910*900 mm
Backrest adjustable	0-75°
Footrest adjustable	0-40°
Height adjustable	450-700mm
Trendeleburg	0-12°
Revere- Trendelenburg	0-12°
Loading capacity	205kg



REMOTE CONTROL

The integrated remote control with clear icons allows nurses to quickly and clearly understand the sections they need to control.



MOTORS

CE mark motors can be replaced according to requirements. Strong stability, great support without any noise. Durable and high safety motor system with compliance to IEC60601 standard, low noise.



DETACHABLE HEADBOARD

The lock structure at the end of the bed is convenient for quick removal in an emergency, allowing medical personnel to be closer to patient's head carry out rescue.



SOFT LINK

Soft link makes platform more smooth to let the patient feel more comfortable when adjust the position.



PP SAFETY RAILING

Tuck-away PP safety railing with damper for slow dropping, special patented design to avoid accidental falling from patient's on bed operation.



CENTRAL LOCKING SYSTEM

Central locking system with dual sided wheels is with excellent stability and convenient for medical staff to manipulate.